





## Sepray CPAP 25/A25

### Comfortable and Effective Solution for OSAHS

Sepray CPAP 25 aims to offer precise pressure to guarantee an effective treatment. Intelligent pressure relief COMF helps to improve the comfort of therapy.



Micomme Sepray CPAP A25 distinguishes between obstructive and central apnea with advanced algorithm. It automatically and continuously adjusts pressure exactly to fit patient's needs to ensure therapeutic effectiveness.



#### User-friendly Design

- 3.5" colour LCD screen and user friendly interface.
- Auto-saving treatment Parameter.

#### Intelligence Core

- With intelligent algorithm, Sepray CPAP A25 reliably differentiates between obstructive and central apnea.
- Sepray CPAP A25 accurately calculates the automatic titration.

#### Safety Care

- Patented 360° anti-countercurrent water chamber.
- Unique anti-dry technology.



#### Data Management

- Smartphone apps and built-in Wi-Fi module make connectivity to Cloud platform in real time.\*
- SpO<sub>2</sub> module can monitor and record blood oxygen saturation and pulse rate.\*



#### Comfortable Treatment

- Powerful pressure relief COMF in all modes.
- Integrated design for main device and humidifier.

Description	Sepray CPAP 25	Sepray CPAP A25
Mode	CPAP	CPAP, Auto CPAP
Pressure range	4 cmH <sub>2</sub> O-20 cmH <sub>2</sub> O	4 cmH <sub>2</sub> O-20 cmH <sub>2</sub> O
Pressure accuracy	± 0.4 cmH <sub>2</sub> O	± 0.4 cmH <sub>2</sub> O
Max. operation pressure	< 30 cmH <sub>2</sub> O	< 30 cmH <sub>2</sub> O
COMF pressure relief	1-3 levels	1-3 levels
Ramp time	0 to 45 min (5-minute increments)	0 to 45 min (5-minute increments)
Humidification	1-5 levels ( 113 to 185°F/45 to 85°C)	1-5 levels ( 113 to 185°F/45 to 85°C)
Humidification settings	Heated Tube (optional)	Heated Tube (optional)
Data storage capacity	> 1 year, Micro SD Card 4G	> 1 year, Micro SD Card 4G
Dimensions (L x W x H)	24.5 x 16 x 10.5 cm	24.5 x 16 x 10.5 cm
Weight	1.82 kg	1.82 kg
Electrical requirements	100 - 240 VAC, 50/60 Hz, 2.0-0.4A	100 - 240 VAC, 50/60 Hz, 2.0-0.4A
Mean sound level	≤ 30 dB	≤ 30 dB
Type of protection against electric shock	Class II Equipment	Class II Equipment
Degree of protection against electric shock	Type BF	Type BF
Degree of protection against ingress of water	IP21	IP21
Relative humidity operation and storage	15% to 95% no condensation	15% to 95% no condensation
Temperature range		
Operation	+5°C to +35°C	+5°C to +35°C
Storage	-25°C to +70°C	-25°C to +70°C



CPAP 25



CPAP A25

Manufacturer: Micomme Medical Technology Development Co., Ltd.

Email: [intl@micomme.com](mailto:intl@micomme.com) Website: [www.micomme.com/en](http://www.micomme.com/en) Tel: +86-4000-2000-33 Fax: +86-0731-89916333

Add: No.8, Lujing Road, High-Tech Zone, Yuelu Changsha, China

NOTE: Specifications subject to change without prior notice.

\* Part of the function is only available for specified countries and areas.

